

For Research Use Only

V5-tag Monoclonal antibody (SV5-P-K)



www.ptglab.com

Catalog Number: v5ab

5 Publications

Basic Information

Catalog Number:
v5ab

Clone No.:
SV5-P-K

Applications:
WB

Host/IsoType:
Mouse/IgG1 kappa

Conjugate:
Unconjugated

Type:
Primary Antibody

Class:
Monoclonal

RRID:
AB_2868500

Purification Method:
Protein A purified

Description

Mouse monoclonal [SV5-P-K] to V5-tag

Specificity/Target

V5-Tag (GKPIPNPLLGLDST)

Recommended Dilution

Western Blot: 1:1,000

Storage

Storage:

Shipped on blue ice. Store at -20°C/-4°F. Stable for 1 year after shipment.

Storage Buffer:


PBS, preservative: 0.02 % Proclin

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Western Blot (WB) analysis of His-eGFP-V5 added to HEK293T cell extract



The anti-V5 antibody [SV5-P-K] is a recombinant, monoclonal mouse IgG1, has very low background and detects even low amounts of V5-tagged proteins. Secondary antibody ChromoTek's Alpaca anti-mouse IgG1 Alexa Fluor 488 Nano-Secondary.