

For Research Use Only

# anti-SNAP/CLIP recombinant VHH, unconjugated



www.ptglab.com

Catalog Number: wt

1 Publications

## Basic Information

### Catalog Number:

wt

### Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

### Host:

Alpaca

### Conjugate:

Unconjugated

### Type :

Nanobody

### Class:

Recombinant

### RRID:

AB\_2827564

### Molecular Weight:

13.1 kDa

## Description

Alpaca anti-SNAP/CLIP-tag VHH, purified recombinant binding protein

## Specificity/Target

SNAP-tag and CLIP-tag

## Affinity ( $K_D$ )

Dissociation constant  $K_D$  of 1 nM

## Storage

### Storage:

Shipped at ambient temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

### Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09% sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com  
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Structure of SNAP/CLIP-tag.