

For Research Use Only

anti-PARP1 recombinant VHH, unconjugated



www.ptglab.com

Catalog Number: xt

Basic Information

Catalog Number:

xt

Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:

Alpaca

Conjugate:

Unconjugated

Type :

Nanobody

Class:

Recombinant

RRID:

AB_2631385

Molecular Weight:

14.36 kDa

Description

Alpaca anti-PARP1 VHH, purified recombinant binding protein

Specificity/Target

Human PARP1. The anti-PARP1 VHH does not bind to human PARP2, 3, or 9
Epitope: within DNA binding domain of PARP1

Storage

Storage:

Shipped at ambient temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.

Selected Validation Data



Structure of PARP1.